


VARIAN INSTRUMENTS AT WORK



ATOMIC ABSORPTION

Multielement Analysis of Rocks and Sediments by Wet Digestion and Atomic Absorption Spectroscopy*

Number AA-47
April 1985

Arthur J. Horowitz and Kent A. Elrick,
U.S. Geological Survey, WRD,
6481-H Peachtree Industrial Blvd.,
Doraville, GEORGIA 30340, U.S.A.

Summary

A relatively rapid, accurate, and precise method for the determination of seventeen elements (Fe, Mn, Al, Ti, Mg, Na, K, Ca, Cu, Zn, Pb, Cr, Ni, Co, Cd, Li, and Sr) in rocks and sediments, from as little as 250 mg of sample is presented. The method employs a wet digestion in Teflon beakers using a combination of HF, HClO₄ and HNO₃ acids. Quantitation is carried out by flame atomic absorption spectrophotometry using mixed salt standards and modifiers. Comparison of data generated by this procedure, with data for USGS Standard Rocks, and NBS sediments, as well as through interlaboratory comparisons on natural sediment samples, indicate that precise and accurate results can be obtained.

* Also published as USGS Open-File Report 85-78.

Introduction

The chemical analysis of rocks and sediments is performed for a variety of environmental and petrological purposes. For example, chemical data are used in the identification of various rock types; also, sediments contain significantly higher concentrations of many metals than are found in the underlying water; as such, they must be analyzed for pollutant contributions to the environment. The importance of these types of data is well established.

Silicate dissolution usually involves either a fusion, with subsequent dissolution of the bead, or a wet digestion employing mineral acids. Fusions tend to be limited to major element analyses due to the relatively high dilution factors involved.

A desire to quantitate a number of trace (heavy) metals dictates the use of a wet digestion procedure. Various wet digestion procedures have been used with rocks and sediments and have been amply described in the literature (1-7). The application of atomic absorption spectroscopy for the analysis of the digestate has also been amply described (1-7). The digestion technique described here represents a modification of the procedures outlined by Johnson and Maxwell (5) and Langhmyr and Paus (6). Table I specifies the upper and

lower concentration limits; samples containing analyte concentrations greater than the upper limit may be analyzed after appropriate dilution.

Summary of Method

Rocks and/or sediment samples are dried, ground and homogenized. An aliquot is digested with a combination of nitric, hydrofluoric, and perchloric acids in Teflon beakers, heated on a hot plate at 200°C. The resulting salts are dissolved in hydrochloric acid and deionized water. The solutions are analyzed by flame atomic absorption spectroscopy after the addition, in certain cases, of appropriate modifiers. Additional interferences are removed and/or compensated for through the use of mixed-salt standards and background correction. Further information about the principles of the method can be found in Walsh (8), Johnson and Maxwell (5) and Pinta (7).

Interferences

Numerous interelement interferences, both positive and negative, exist for this procedure and have been amply documented elsewhere (1,5,7,8). Interferences are eliminated and/or compensated for through the removal of silica by the digestion procedure, dilution, the addition of cesium chloride (CSCl), the use of mixed-salt standards, and background correction.

Apparatus

A Varian Model AA-975¹¹ double beam atomic absorption spectrophotometer with microprocessor control and digital display, used in conjunction with a Varian Model PSC-55¹¹ autosampler was employed in this study. Instrumental parameters are listed in Table II.

Teflon beakers, 100 mL capacity, thick wall, capable of withstanding temperatures up to 260°C.

Hot plate, electric or gas, capable of at least 250°C.

Perchloric acid hood, with appropriate washdown facility and gas or electric outlets.

Reagents

Aluminium standard solution, 1.00 mL = 1.00 mg Al: Dissolve 1.000 g aluminium metal in 20 mL HCl (sp gr 1.19) with a trace of mercury salt to catalyze the reaction, and dilute to 1000 mL with demineralized water.

Cesium chloride solution, 4 g/L: Dissolve 4 g CsCl of at least 99.999% purity in demineralized water and dilute to 1 L.

Hydrochloric acid, concentrated, (sp gr 1.19).

Hydrochloric acid, dilute (1 + 1): Add 250 mL concentrated hydrochloric acid (sp gr 1.19) to 250 mL demineralized water. Store in a plastic bottle.

Hydrochloric acid, dilute (2 + 98): Add 10 mL concentrated hydrochloric acid (sp gr 1.19) to 490 mL demineralized water. Store in a plastic bottle.

Hydrofluoric acid, concentrated (48-51%), (sp gr 1.17).

Iron standard solution, 1.00 mL = 1.00 mg Fe: Dissolved 1.000 g iron metal in 20 mL HCl (1 + 1) and dilute to 1000 mL with demineralized water.

Mixed salt standard stock solution I (minors):

Dissolve by appropriate means the following compounds or elements: cadmium metal (0.200 g), chromium metal (0.800 g), cobalt metal (1.200 g), copper metal (0.800 g), lead metal (2.000 g), lithium carbonate (2.130 g), manganese metal (2.000 g), nickel metal (1.200 g), strontium carbonate (1.685 g), and zinc metal (0.320 g); add 20 mL HCl (sp gr 1.19) and dilute to 1000 mL with demineralized water. This solution will contain the following concentrations:

cadmium (200 mg/L), chromium (800 mg/L), cobalt (1200 mg/L), copper (800 mg/L), lead (2000 mg/L), lithium (400 mg/L), manganese (2000 mg/L), nickel (1200 mg/L), strontium (1000 mg/L), and zinc (320 mg/L). Store in a plastic or Teflon bottle.

Mixed salt standard stock solution IA (minors):

Take 100 mL of mixed salt standard stock solution I, add 20 mL HCl (sp gr 1.19), and dilute to 1000 mL. This solution will contain the following concentrations:

cadmium (20 mg/L), chromium (80 mg/L), cobalt (120 mg/L), copper (80 mg/L), lead (200 mg/L), lithium (40 mg/L), manganese (200 mg/L), nickel (120 mg/L), strontium (100 mg/L), and zinc (32 mg/L). Store in a plastic or Teflon bottle. Solution is stable for 3 months.

Mixed salt standard solution II (majors):

Dissolve by appropriate means the following compounds or elements: aluminium metal (1.500 g), calcium carbonate (1.249 g), iron metal (1.000 g), magnesium metal (0.200 g), manganese metal (0.040 g), potassium chloride (0.668 g), sodium chloride (0.636 g), and ammonium titanyl oxalate (1.227 g); add 20 mL HCl (sp gr 1.19) and dilute to 1000 mL with demineralized water.

This solution will contain the following concentrations: aluminium (1500 mg/L), calcium (500 mg/L), iron (1000 mg/L), magnesium (200 mg/L) manganese (40 mg/L), potassium (350 mg/L), sodium (250 mg/L) and titanium (200 mg/L). Store in a plastic or Teflon bottle.

Working standard solution I:

Take respectively, a 10 mL, 5 mL and 1 mL aliquot of mixed-salt standard stock solution IA, add to each 4 mL HCl (sp gr 1.19), 20 mL of mixed-salt standard stock solution II, and dilute to 200 mL in volumetric glassware with demineralized water. Store in plastic or Teflon bottles. Prepare fresh for each analysis. Concentrations are as follows:

	Standard 3 (mg/L)	Standard 2 (mg/L)	Standard 3 (mg/L)
Volume (mL)	10	5	1
Cd	1	0.5	0.1
Cr	4	2	0.4
Co	6	3	0.6
Cu	4	2	0.4
Pb	10	5	1
Li	2	1	0.2
Ni	6	3	0.6
Sr	5	2.5	0.5
Zn	1.6	0.8	0.16

Working standard solution II:

Take respectively, a 10 mL, 6 mL, and 2 mL aliquot of mixed-salt standard stock solution II, add 2 mL HCl (sp gr 1.19) and 10 mL of the CsCl solution, and dilute to 100 mL in volumetric glassware with demineralized water. Store in plastic or Teflon bottles. Prepare fresh for each analysis.

Concentrations are as follows:

	Standard 4 (mg/L)	Standard 5 (mg/L)	Standard 6 (mg/L)
Volume (mL)	10	6	2
Al	150	90	30
Fe	100	60	20
Mg	20	12	4
Mn	4	2.4	0.8

Working standard solution III:

Take a 10 mL aliquot of standards 4, 5 and 6, add 2 mL HCl (sp gr 1.19), and 10 mL of the CsCl solution and dilute to 100 mL in volumetric glassware with demineralized water. Store in plastic or Teflon bottles. Prepare fresh for each analysis. Concentrations are as follows:

	Standard 7 (mg/L)	Standard 8 (mg/L)	Standard 9 (mg/L)
Volume (mL)	10 (Standard 4)	10 (Standard 5)	10 (Standard 6)
Ca	5.0	3.0	1.0
K	3.5	2.1	0.7
Na	2.5	1.5	0.5

Nitric acid, concentrated (sp gr 1.41).

Perchloric acid, concentrated (70-72%) (sp gr 1.67).

Sodium standard solution, 1.00 mL = 1.00 mg Na:
Dissolve 2.542 g NaCl in demineralized water, add 20 mL HCl (sp gr 1.19), and dilute to 1000 mL with demineralized water.

Titanium standard solution, 1.00 mL = 1.00 mg Ti:
Dissolve 6.135 g of ammonium titanyl oxalate in demineralized water, and dilute to 1000 mL with demineralized water.

Titanium working standard solutions:

Take respectively, a 2 mL, 1 mL and 0.5 mL aliquot of the titanium standard solution, add to each 10 mL of the aluminium standard solution, 5 mL of the iron standard solution, 3.5 mL of the sodium standard solution, 10 mL of the CsCl solution and 2 mL HCl (sp gr 1.19) and dilute to 100 mL in volumetric glassware with demineralized water. The standards contain, respectively, 20, 10 and 5 mg/L titanium.

Procedure

Immediately before each use, clean all glassware by rinsing, first with dilute HCl (1 + 1), and then with demineralized water. Dry the sample by an appropriate procedure such as freeze-drying, or in an oven at 105°C. If the sample is greater than 100 g, split it down to less than 100 g by use of a non-metallic sample splitter (riffle sampler) or by coning and quartering. Grind the sample with a mixer mill or an agate mortar and pestle until all material is finer than 100 mesh.

Weigh and transfer 0.5000 g of finely ground sample to a 100 mL Teflon beaker; weigh out appropriate standard materials as well, and use several empty beakers for blanks^[2].

Place the hot plate in a perchloric acid hood, turn on the hood and hot plate, and adjust the hot plate to produce a surface temperature of 200°C. To each beaker, add 6 mL HNO₃, (sp gr 1.41) and place it on the hot plate for approximately 30 minutes^[3].

Remove the beakers from the hot plate and wait 5 minutes. Add 6 mL HF (sp gr 1.17) and 2 mL HClO₄ (sp gr 1.67), and return the beakers to the hot plate. Continue heating until the evolution of white perchloric fumes and the solutions have reached incipient dryness; however, do not bake the residues. Remove the beakers from the hot plate, wait 5 minutes, and repeat the process again. Remove the beakers from the hot plate, wait 5 minutes, and add 2 mL HClO₄ (sp gr 1.67) and return the beakers to the hot plate. Continue heating until the evolution of white perchloric fumes and the solution reaches incipient dryness; however, do not bake the residues. Remove the beakers from the hot plate, lower the hot plate temperature to 100°C, and add 2 mL dilute HCl (1+1) and swirl the beaker; add 10 mL demineralized water and return to the hot plate until the residues dissolve.

Cool the beakers, and pour each solution into a 50 mL volumetric flask. Rinse the beaker several times with demineralized water and bring to the mark with demineralized water^[4]. Pour the solution into an acid-rinsed plastic bottle for storage. This solution represents a dilution factor of 100X. Remove a 5 mL aliquot from the 100X solution, add 1 mL HCl (sp gr 1.19), and 5 mL CsCl solution^[5], place in a 50 mL volumetric flask, and bring to the mark with demineralized water. Pour the solution into an acid-rinsed plastic bottle for storage. This solution represents a dilution factor of 1000X. Finally, remove a 5 mL aliquot from the 1000X solution, add 1 mL HCl (sp gr 1.19), and 5 mL CsCl solution^[5], place in a 50 mL volumetric flask, and bring to the mark with demineralized water. Pour the solution into an acid-rinsed plastic bottle for storage. This solution represents a dilution factor of 10,000X.

Set up the atomic absorption spectrophotometer according to the specifications outlined in Table II, and analyze the 100X solutions for Cd, Cr, Co, Cu, Pb, Ni, and Zn using standards 1, 2 and 3. Dilute samples further if required. Transfer 5 mL aliquots of each sample and standard to an appropriate container, add 5 mL dilute HCl (2 + 98), and 1 mL CsCl solution, and analyze the solutions for Li and Sr using the conditions listed in Table II^[6].

Set up the atomic absorption spectrophotometer as outlined in Table II, and analyze the 1000X solutions for Fe, Mn, Mg and Al using standards 4, 5 and 6. Dilute further if required. Also analyze the 1000X solutions for Ti using the Ti working standards^[7].

Set up the atomic absorption spectrophotometer as outlined in Table II, and analyze the 10,000X solutions for Ca, K and Na using standards 7, 8 and 9. Dilute samples further if required.

Calculations and Reporting Limits

Determine the concentration of each constituent in the 100X solutions (Cd, Cr, Co, Cu, Pb, Ni and Zn) while aspirating each sample and record the results (average of 3 readings once the system has stabilized). The actual concentration of each constituent in the sample can be obtained by multiplying the concentration in each sample solution by 100, if no further dilutions are made.

Determine the concentration of each constituent in the 200X solutions (Li and Sr) while aspirating each sample and record the results (average of 3 readings once the system has stabilized). The actual concentration of each constituent in the sample can be obtained by multiplying the concentration in each sample solution by 200, if no further dilutions are made.

Determine the concentration of each constituent in the 1000X solutions (Fe, Mn, Mg, Al and Ti) while aspirating each sample and record the results (average of 3 readings once the system has stabilized). The actual concentration of each constituent in the sample can be obtained by multiplying the concentration in each sample solution by 1000, if no further dilutions are made.

Determine the concentration of each constituent in the 10,000X solutions (Na, K, and Ca) while aspirating each sample and record the results (average of 3 readings once the system has stabilized). The actual concentration of each constituent in the sample can be obtained by multiplying the concentration in each sample solution by 10,000, if no further dilutions are made.

The reporting limits for each major constituent are as follows: aluminium (nearest 1000 mg/kg), calcium (nearest 1000 mg/kg), iron (nearest 1000 mg/kg), magnesium (nearest 1000 mg/kg), manganese (nearest 100 mg/kg), potassium (nearest 1000 mg/kg), sodium (nearest 1000 mg/kg), and titanium (nearest 1000 mg/kg). As 10,000 mg/kg equals 1%, Al, Ca, Fe, Mg, K, Na and Ti should be reported to the nearest tenth of a percent; Mn should be reported to the nearest hundredth of a percent. The reporting limits for each minor constituent are as follows: cadmium (nearest 0.1 mg/kg to 10, above 10, nearest mg/kg), chromium (nearest 1 mg/kg), cobalt (nearest 1 mg/kg), copper (nearest 1 mg/kg), lead (nearest 1 mg/kg), lithium (nearest 1 mg/kg), nickel (nearest 1 mg/kg), strontium (nearest 1 mg/kg) and zinc (nearest 1 mg/kg).

Precision and Accuracy

The precision and accuracy of this method was determined by replicate analyses (actual separate digestions and subsequent quantitation) on 2 National Bureau of Standards Standard Reference Materials, and 6 U.S. Geological Survey Rock Standards. The results are presented in Table III. As can be seen from the data, the method is capable of generating both precise and accurate analytical results.

Discussion and Conclusions

In order to further evaluate the precision and accuracy of this method, 17 natural freshwater and marine sediment samples were dried, digested and analyzed. The samples came from different geological settings and water bodies (Appalachicola River, Florida; Patuxent River, Maryland; Doane Lake, Swan Island and Columbia Slough, Oregon; Mississippi River, Louisiana; Ned Wilson Lake, Colorado; Yaharra and Nemadji Rivers, Wisconsin; Lake Bruin, Louisiana; and George's Bank from the North Atlantic Outer Continental Shelf). The results for the major element analyses were compared with the results obtained from a fusion procedure with subsequent atomic absorption quantitation (Table IV). The results from the two sets of analyses are remarkably consistent.

As an additional check for the minor elements, 6 of the samples, already dried and ground, were sent to another laboratory for dissolution and subsequent quantitation. The results of the two sets of analyses are presented in Table V. Quantitation by the outside laboratory for all elements but Ti, came from an HF/HClO₄/HNO₃ digestion with subsequent quantitation by flame or furnace AA. Ti was determined by ICP-AES following a borate fusion.

As can be seen from all the comparative data on standards and samples, analytical precision and accuracy, as well as comparability, is quite good (Tables III, IV and V). All these results indicate that very precise and accurate analyses can be obtained on rocks and sediments by using flame atomic absorption spectroscopy.

References

1. Abbey, S., Canadian Mineralogist, **8**, 347-353 (1965).
2. Agemian, H. and Chau, A., Analytica Chimica Acta, **80**, 61-66 (1975).
3. Bernas, B., Analytical Chemistry, **40**, 1682-1686 (1968).
4. Dolezal, J., Povondra, P., and Sulcek, Z., Decomposition Techniques in Inorganic Analysis, Elsevier Publishing Co., Inc., New York, NY., 11-89 (1966).
5. Johnson, W., and Maxwell, J., Rock and Mineral Analysis, 2nd ed., John Wiley and Sons, New York, NY., 237-325 (1981).
6. Langmyhr, F., and Paus, P., Analytica Chimica Acta, **43**, 397-408 (1968).
7. Pinta, M., Modern Methods for Trace Element Analysis, Ann Arbor Science Publishers, Ann Arbor, MI., 208-219 (1982).
8. Walsh, J., Analyst, **102**, 972-976 (1977).

Footnotes

- [1] The use of brand names in this report is for identification purposes only, and does not constitute an endorsement by the U.S. Geological Survey.
- [2] This procedure can be used with sample weights of between 0.2500 and 1.0000 g, with appropriate adjustments to the final solution volumes and acid strengths. Larger sample weights may be used, but will almost certainly require a triple digestion with HF and HClO₄.
- [3] This step is designed to oxidize organic matter in the sample. It is imperative that this step be carried out prior to the addition of perchloric acid, otherwise a violent explosion could occur.
- [4] If a sample contained a large amount of organic matter, it is not unusual to have black 'flecks' in the final solution; these can be ignored if allowed to settle prior to aspiration into the AAS.
- [5] The CsCl acts as an ionization suppressant.
- [6] The additional dilution is required to eliminate interferences due to density differences (1), and the CsCl acts as an ionization suppressant.
- [7] Titanium determinations by atomic absorption are subject to severe interferences and sensitivity is heavily dependent on flame stoichiometry (6,8). Adjust the nitrous oxide flame so that it is nearly luminous (increase the fuel flow until the reducing red cone turns orange-yellow, then reduce the fuel flow until the flame just becomes red again). The CsCl acts as an ionization suppressant.

TABLE I

Upper and Lower Concentration Limits of the Method

Constituent	Lower Limit		Upper Limit	
	Sample Solution (mg/kg)	Solution (mg/kg)	Sample Solution (mg/kg)	Solution (mg/kg)
Aluminium	20,000	20	150,000	150
Cadmium	0.5	0.005	10	0.1
Calcium	1,000	0.1	50,000	5
Chromium	3	0.03	400	4
Cobalt	3	0.03	600	6
Copper	1	0.01	400	4
Iron	5,000	5	100,000	100
Lead	3	0.03	1,000	10
Lithium	2	0.02	200	2
Magnesium	1,000	1	20,000	20
Manganese	100	0.1	4,000	4
Nickel	3	0.03	600	6
Potassium	1,000	0.1	35,000	3.5
Sodium	1,000	0.1	25,000	2.5
Strontium	5	0.025	1,000	5
Titanium	1,000	1	20,000	20
Zinc	1	0.01	160	1.6

TABLE II

Instrumental Settings Used for the Method

Parameters	Al	Ca	Fe	Mg	Mn	K	Na	Ti	Cu	Zn	Pb	Ni	Co	Cr	Cd	Li	Sr
Wavelength (nm)	309.3	422.7	372.0	202.6	279.5	766.5	589.0	364.3	324.7	213.9	217.0	232.0	240.7	357.9	228.8	670.8	460.7
Slit (nm)	0.5	0.5	0.2	1.0	0.2	1.0	1.0	0.5	0.5	1.0	1.0	0.2	0.2	0.2	0.5	1.0	0.5
Lamp Current (mA)	10	4	5	4	5	5	5	20	4	5	5	4	7	7	4	5	10
Flame Type (a,b)	N-A	N-A	A-A	N-A	A-A	A-A	A-A	N-A	A-A	A-A	A-A	A-A	A-A	N-A	A-A	A-A	N-A
Oxidant (L/min) ^c	12.8	12.0	12.5	11.0	14.0	14.0	16.2	11.0	11.0	11.0	11.5	11.0	11.0	11.0	11.0	11.5	12.0
Fuel (L/min) ^c	7.00	6.50	2.00	7.00	2.10	3.10	2.20	7.20	2.20	2.20	2.20	2.20	2.20	6.50	2.20	2.20	5.50
Integration Time (s)	5	3	3	5	3	3	3	5	3	3	3	3	3	5	3	3	5
Background Corr.	off	off	off	off	off	off	off	off	on	on	on	on	on	off	on	off	off

*All determinations carried out with a fixed-rate nebulizer, approximate uptake of 5-6 ml/minute.

- a. N-A = nitrous oxide-acetylene
b. A-A = air-acetylene
c. Readings on automatic gas control

TABLE III

Comparison of Reported Concentration with those Found by Method

Element	No. of Det's	WBS River Sediment	WBS Zatuarina Sediment	W-2	C-2	BCE-1	MAG-1	SCR-1	SCO-1
Cu	rpt ¹ fnd ²	109±19 106±2	18±3 18±1	106±5 105±2	10±2 10±2	16±3 16±1	33±3 31±1	65±3 62±2	30±3 29±2
Zn	rpt fnd	1720±170 1 3 4 ± 3	138±6 134±3	80±3 80±3	88±4 85±1	127±4 130±2	132±7 130±3	88±6 85±1	110±6 106±2
Pb	rpt fnd	714±28 708±7	28.2±1.8 28±3	14±2 14±1	29±2 29±2	14±2 13±1	25±2 25±1	40±3 41±1	28±3 29±2
Ni	rpt fnd	45.8±2.9 43±2	32±3 30±1	67±3 66±2	3±1 3±1	10±2 9±1	52±3 51±1	31±3 31±2	28±3 27±1
Co	rpt fnd	10.1±0.6 9±1	10.5±1.3 12±1	43±2 44±2	5±1 6±1	39±2 40±2	19±2 22±1	11±2 13±1	11±2 12±1
Cd	rpt fnd	10.2±1.5 10±0.04	0.36±0.07 0.4±0.08	<1 0.5	<1 0.5	<1 0.5	<1 0.5	<1 0.5	<1 0.5
Li	rpt fnd	— —	49 47±1	10±1 8±1	34±2 34±1	14±2 14±1	75±1 73±1	1 3 ± 5 132±5	44±2 45±2
Cr	rpt fnd	— —	76±3 75±2	92±4 91±3	8±1 8±1	14±2 13±1	103±3 102±3	31±2 31±1	65±2 65±3
Sr	rpt fnd	— —	— —	192±3 208±14	480 475±14	330 341±8	150 144±10	430 430±13	175 163±6
Fe	rpt fnd	11.3±1.2 10.8±0.2	3.35±0.1 3.21±0.05	7.6±0.2 7.5±0.1	1.9 1.8±0	9.4 9.4±0.2	4.8 4.6±0.05	2.2 2.1±0.07	3.6 3.6±0.1
Mn	rpt fnd	0.08±0.01 0.07±0.005	0.04±0.002 0.04±0	0.13±0.01 0.13±0	0.02 0.02±0.005	0.14 0.14±0	0.08 0.08±0.004	0.02 0.02±0	0.04 0.04±0
Mg	rpt fnd	0.074±0.02 0.7±0	1.09±0.08 1.1±0.05	3.8±0.1 3.8±0.15	0.45 0.44k±0.01	2.1 2.0±0.04	1.8 1.7±0.04	2.7 2.6±0.05	1.6 1.5±0.05
Al	rpt fnd	2.26±0.04 2.4±0.04	6.2±50.2 6.3±0.15	8.2±0.1 8.2±0.11	8.2 8.3±0.11	7.3 7.2±0.18	8.7 8.7±0.25	3.5 3.4±0.08	7.2 7.2±0.15
Ti	rpt fnd	— —	0.51 0.46±0.03	0.67±0.02 0.68±0.02	0.3 0.3±0.02	1.4 1.4±0.06	0.4 0.4±0.002	0.2 0.2±0.01	0.4 0.4±0.03
Ca	rpt fnd	2.9 2.8±0.09	0.83±0.03 0.74±0.05	7.8±0.1 7.7±0.08	1.4 1.3±0.04	5.0 4.9±0.08	1.0 0.9±0	6.4 6.3±0.08	1.9 1.8±0.04
Na	rpt fnd	0.54±0.01 0.56±0.09	2.0 2.0±0.09	1.6±0.01 1.6±0.04	3.0 2.9±0.1	2.5 2.5±0.05	2.8 2.8±0.05	2.2 2.2±0.09	0.7 0.7±0.07
K	rpt fnd	1.26±0.05 1.26±0.05	1.4 1.5±0.05	0.5±0.01 0.6±0.05	3.7 3.6±0.07	1.4 1.4±0	3.0 2.9±0.05	1.3 1.3±0.04	2.2 2.2±0

¹rpt:reported

²fnd:found

TABLE IV
Comparison of Major Element Concentrations Determined by Fusion and by This Method

	Fe ¹		Mn ¹		Mg ¹		Al ¹		Ca ¹		Na ¹		K ¹	
	F ²	W	D ³	F	WD	F	WD	F	WD	F	WD	F	WD	F
Appalachicola River, R.M. 94	1.1	1.1	0.04	0.04	0.06	0.07	1.9	2.1	0.26	0.25	0.4	0.5	1.0	1.0
Patuxent River at Hog Point	2.8	3.0	0.04	0.05	0.75	0.74	5.1	4.9	0.5	0.4	1.3	1.4	1.6	1.5
Doane Lake Outlet	6.0	6.1	0.08	0.09	1.0	1.0	8.1	8.0	2.3	2.2	1.9	1.8	1.3	1.2
Mississippi River at Venice	3.7	3.8	0.09	0.09	1.0	1.1	7.6	7.6	1.0	1.0	0.8	0.9	2.1	2.2
Ned Wilson Lake	2.8	2.8	0.04	0.03	1.3	1.3	7.2	7.2	1.0	1.0	0.8	0.9	1.9	1.9
Yaharra River	1.2	1.3	0.07	0.06	2.2	2.3	3.3	3.3	4.7	4.6	0.3	0.5	3.0	3.2
Appalachicola River, R.M. 86	0.05	0.10	0.02	0.01	0.01	0.02	0.3	0.4	0.1	0.1	0.5	0.4	0.3	0.3
Appalachicola River, Sandbar	0.25	0.29	0.02	0.02	0.02	0.02	0.4	0.4	0.1	0.1	0.3	0.3	0.3	0.4
Mississippi River at Tarbert	0.48	0.48	0.01	0.01	0.09	0.09	1.5	1.6	0.4	0.4	1.0	0.9	1.0	1.1
George's Bank M8 5-4	0.48	0.43	0.01	0.02	0.07	0.07	0.2	0.3	0.2	0.2	0.6	0.5	0.3	0.2
Patuxent River at Pt. Patience	2.2	2.3	0.05	0.06	0.5	0.6	3.9	4.0	0.6	0.5	1.2	1.0	1.2	1.2
Patuxent River at St. Leonard's Ck.	3.8	4.1	0.08	0.10	0.8	0.9	7.2	7.2	0.6	0.6	1.2	1.1	1.3	1.5
Columbia Slough	3.5	3.7	0.07	0.06	1.0	1.0	7.9	8.1	2.6	2.5	2.1	2.2	1.6	1.5
Swan Island	4.9	4.8	0.09	0.09	1.2	1.2	8.5	8.5	2.6	2.5	2.0	2.1	1.3	1.4
Nemadji River	2.4	2.4	0.07	0.06	1.1	1.1	4.3	4.3	1.8	1.7	0.9	1.0	2.5	2.5
George's Bank M13A	2.4	2.5	0.05	0.04	0.7	0.8	4.8	4.9	1.0	0.9	1.8	2.0	2.0	1.8
Lake Bruin	3.3	3.2	0.10	0.09	0.7	0.8	6.6	6.5	0.8	0.8	1.0	0.9	2.0	1.8

x¹ — all concentrations in weight percent (1% = 10,000 mg/kg).
 F² — fusion with lithium metaborate/tetraborate followed by AAS quantitation.
 WD³ — this method.

TABLE V
Interlaboratory Comparison of Minor Element Data in Selected Sediment
(concentration in mg/kg)

Sample		Cu	Zn	Pb	Cr	Cd	Ni	Co	Li	Sr	Ti
Appalachicola River RM 94	a ¹	4	23	12	23	<0.5	4	10	9	43	4500
	b ²	6	23	15	20	0.03	6	10	8	ND ³	4600
Patuxent River at Hog Point	a ¹	20	111	22	56	0.8	26	32	40	60	4700
	b ²	20	110	25	59	0.5	27	33	41	50	4900
Doane Lake Outlet	a ¹	33	127	29	58	<0.5	34	25	16	300	6600
	b ²	35	130	32	56	0.14	29	28	17	290	6600
Mississippi River at Venice	a ¹	28	120	32	76	0.9	37	16	38	100	4800
	b ²	30	120	32	77	0.5	36	16	39	90	4900
Ned Wilson Lake	a ¹	28	113	38	90	0.9	66	20	36	210	4900
	b ²	30	110	37	90	0.4	55	20	37	190	4700
Yaharra River	a ¹	13	28	23	28	0.5	11	18	6	107	1900
	b ²	12	27	22	30	0.1	13	16	7	94	1900

a¹ — this method
 b² — USGS Branch of Analytical Services, Reston, VA, determined on a sample digested with HF/HClO₄/HNO₃
 Cu — digestate extracted with butyl acetate, further extracted with diethyldithiocarbamate in chloroform, dried, brought up in HCl, graphite furnace
 Zn — see Cu, quantitation by flame AAS
 Pb — see Cu
 Cd — see Cu
 Cr — digestate diluted 1:10, quantitation by graphite furnace
 Ni — see Cr
 Co — see Cr
 Li — digestate diluted, quantitation by flame AAS
 Sr — see Li
 Ti — borate fusion, quantitation by ICP
 ND³ — not determined