

Varian Laboratory Services

Reliable scientific support to help secure your world

For more than 12 years, Varian's state-of-the-art drug testing laboratory has provided confidential and convenient drug-testing services that help detect and deter substance abuse. Our full range of testing options, backed by Varian's commitment to customer care and quality performance, will give you prompt, accurate, and legally defensible results.

Services

Confirmation testing

Varian onsite DAT products provide only a preliminary analytical test result. The preferred confirmation method, Gas Chromatography/Mass Spectrometry (GC/MS), should be used in order to obtain a confirmed analytical result. Varian Lab Services uses either GC/MS or GC/MS/MS to deliver confirmation results within 24–48 hours. Results can be delivered in a variety of ways, including secure fax, phone, and/or email, depending on customer preference.

Virtual Laboratory (V•Lab™) Services **IN DEVELOPMENT**

Rapid laboratory-certified test results from any location

When your onsite urine tests require laboratory-certified results, look to V•Lab. Scan the Varian testing device you prefer into Varian's exclusive TesTrak™ ID data management program. A simple click automatically transmits the secure image to our laboratory. One of our scientists will interpret negative test results in less than four hours.* Confirmation testing results will be returned within two working days. V•Lab meets College of American Pathologists (CAP) requirements for forensic urine drug testing.

Drug screening

- **Urine testing** – reliable and accurate detection of commonly abused drugs.
- **Oral fluid testing** – identification of "under the influence" drug use. Choose Varian's OnTrak OraTube™ and discover a unique new approach that combines the convenience of onsite oral fluid collection with laboratory-certified results and doesn't require subjective sample dilution prior to shipping.
- **Hair testing** – non-invasive and resistant to adulteration, hair testing provides a historical "picture" of drug abuse – up to 90 days prior to testing date.

*A V•Lab specimen received by Varian Laboratory Services before 12 noon CST is guaranteed a turnaround time of four hours or less. A V•Lab specimen received after 12:00 noon CST is guaranteed a turnaround time of less than 24 hours.



VARIAN

**Rapid drug test solutions
that secure your world™**



Additional Services

Substance Abuse and Mental Health Services Administration (SAMHSA)-certified test results are available upon request.

Adulteration testing

As a standard validity check, Varian measures urine creatinine concentration on all specimens to identify diluted or adulterated samples. Testing for specific gravity, pH, and nitrite are available upon request, as is specialized adulterant detection using reagents for chromate, glutaraldehyde, and general oxidants.

Expanded panels

Varian's drug immunoassays are available as expanded panels, allowing a larger aggregation of discrete drug assays while maintaining the immunoassay's high sensitivity.

TOXI•LAB®

Varian's TOXI•LAB system uses modified thin layer chromatography to detect and identify 500 different drugs and metabolites from just one sample. It is ideal for cases that require both comprehensive screening and superior discriminatory characteristics.

Certifications

Varian Laboratory Services maintains all required certifications, state licenses, and accreditations, including:

- Clinical Laboratory Licensure (State of New York)
- Clinical Laboratory Licensure (State of Florida)
- Clinical Laboratory Licensure (State of Illinois)
- Clinical Laboratory Improvement Amendment of 1988 (CLIA) certification for interstate commerce
- College of American Pathologists (CAP) forensic urine drug testing
- Ongoing participation in blind urine proficiency testing performed by CAP and by New York State and Pennsylvania's Departments of Health. Consistently receive a qualitative rating of 100%
- Participation in SAMHSA trial proficiency program for alternative specimens
- Prior reference laboratory for the Society of Hair Testing Proficiency Specimens



VARIAN

25200 Commercentre Drive
Lake Forest, CA 92630-8810

To learn more about Varian Laboratory Services, visit us on the Web at www.varian-onsite.com. Or, call us toll-free at **800•737•9667 x5**.

Inspiring Excellence